

"EINE SICHERE PARTNERSCHAFT  
FÜR IHRE GESUNDHEIT"

**mediland**  
products & more

products & more  
**mediland**



### Thorax-Drain

- aus PVC
- Länge 50cm
- 5 seitliche Augen
- graduiert bei 10, 20, 25, 30cm
- x-ray
- zentrale Öffnung
- steril
- VE 25 Stück
- UK 100 Stück

<b>ML-790016</b>	Ch.16
<b>ML-790018</b>	Ch.18
<b>ML-790022</b>	Ch.22
<b>ML-790024</b>	Ch.24
<b>ML-790026</b>	Ch.26
<b>ML-790028</b>	Ch.28
<b>ML-790030</b>	Ch.30
<b>ML-790032</b>	Ch.32
<b>ML-790034</b>	Ch.34
<b>ML-790036</b>	Ch.36

### Thoracic catheter

- of PVC
- length 50cm
- 5 side eyes
- market at 15, 20, 25, 30cm
- x-ray
- central opening
- sterile
- unit 25 pieces
- box 100 pieces

mediland GmbH • Gewerbepark Rudersberg • Ziegeleistrasse 4 • D-73635 Rudersberg (Germany)  
Tel.: +49 (0) 7183 30 78 01 0 • Fax: +49 (0) 7183 30 78 01 11 • [www.mediland.de](http://www.mediland.de) • [info@mediland.de](mailto:info@mediland.de)

"EINE SICHERE PARTNERSCHAFT  
FÜR IHRE GESUNDHEIT"

**mediland**  
products & more

products & more  
**mediland**



### Thorax-Drain mit Trocar

- 2 seitliche Augen
- Länge 40cm
- steril
- VE 1 Stück
- UK 1 Stück

### Thoracic catheter with trocar

- with 2 side eyes
- length 40cm
- sterile
- unit 1 piece
- box 1 piece

- ML-797016** Ch.16
- ML-797018** Ch.18
- ML-797022** Ch.22
- ML-797025** Ch.25
- ML-797028** Ch.28
- ML-797032** Ch.32
- ML-797036** Ch.36

mediland GmbH • Gewerbepark Rudersberg • Ziegeleistrasse 4 • D-73635 Rudersberg (Germany)  
Tel.: +49 (0) 7183 30 78 01 0 • Fax: +49 (0) 7183 30 78 01 11 • [www.mediland.de](http://www.mediland.de) • [info@mediland.de](mailto:info@mediland.de)