

"EINE SICHERE PARTNERSCHAFT
FÜR IHRE GESUNDHEIT"

mediland
products & more

products & more
mediland



Tiemann-Katheter soft flexibel

- Länge 40cm
- gebogene, flexible Spitze
- 2 seitliche Augen
- Augen poliert
- steril
- VE 10 Stück
- UK 5x10 Stück

Tiemann catheter soft flexible

- length 40cm
- bent tip flexible
- 2 eyes
- polished eyes
- sterile
- unit 10 pieces
- box 5x10 pieces

ML-120FS12	Ch.12
ML-120FS14	Ch.14
ML-120FS16	Ch.16
ML-120FS18	Ch.18
ML-120FS20	Ch.20
ML-120FS22	Ch.22
ML-120FS24	Ch.24



mediland GmbH • Im Fuchshau 22 • D-73635 Rudersberg (Germany)

Tel.: +49 (0) 7183 30 78 01 0 • Fax: +49 (0) 7183 30 78 01 11 • www.mediland.de • info@mediland.de

"EINE SICHERE PARTNERSCHAFT
FÜR IHRE GESUNDHEIT"

mediland
products & more

products & more
mediland



Tiemann-Katheter

- Länge 40cm
- geschlossene Spitze
- 2 seitliche Augen
- steril
- VE 100 Stück
- UK 500 Stück

Tiemann catheter

- length 40cm
- Tiemann tip
- 2 eyes
- sterile
- unit 100 pieces
- box 500 pieces

ML-120006	Ch.6
ML-120008	Ch.8
ML-120010	Ch.10
ML-120012	Ch.12
ML-120014	Ch.14
ML-120016	Ch.16
ML-120018	Ch.18
ML-120020	Ch.20

mediland GmbH • Im Fuchshau 22 • D-73635 Rudersberg (Germany)

Tel.: +49 (0) 7183 30 78 01 0 • Fax: +49 (0) 7183 30 78 01 11 • www.mediland.de • info@mediland.de