


"EINE SICHERE PARTNERSCHAFT  
FÜR IHRE GESUNDHEIT"

**mediland**  
products & more

products & more  
**mediland**



### Ureterokutaneostomie Katheter

- aus weichem Polyurethan
- transparent
- zentrale Öffnung
- cm-Graduierung
- Drainageöffnungen alle 10mm
- Selbstarretierung
- Fixierscheibe
- steril
- VE 2 Stück
- UK 50 Stück

### Catheter for ureter skin fistula

- made of soft Polyurethane
- transparent
- central open tip
- graduated every cm
- perforated every 10mm
- fixation against
- dislocation
- antegrad as well as
- retrograd
- sterile
- unit 2 pieces
- box 50 pieces

**ML-138723** 26cm 3,5mm

**ML-138724** 26cm 4,2mm

**ML-138733** 32cm 3,5mm

**ML-138734** 32cm 4,2mm