


"EINE SICHERE PARTNERSCHAFT
FÜR IHRE GESUNDHEIT"

mediland
products & more

products & more
mediland



Uretersplint-Set „Mono-J“

- Polyurethan (PUR)
- super x-ray
- Länge 90cm
- Ch.4,8 - Ch.9,0
- geschlossene Spitze und zentrale Spitze
- im Spitzenbereich perforiert
- Führungsdraht 100cm - 150cm
- Führungsdraht beschichtet flexible Spitze
- mit Fixierklemmen
- je 1 Katheter mit roter- und grüner Markierung im Set
- steril
- VE 1 Set
- UK 25 Sets

Ureteral splint (set) „Mono-J“

- Polyurethane (PUR)
- super x-ray
- length 90cm
- Ch.4,8 - Ch.9,0
- cylindrical oder central opening
- perforated in the tip area
- guide-wire 100cm - 150cm
- guide-wire coated flexible tip
- with fastening clips
- 1 catheter red and green market in the set
- sterile
- unit 1 set
- box 25 sets

Größe/size

Führungsdraht / guide-wire

**zylindrisch
cylindrical**

**zentral offen
central opening**

Ch.4,8 (0.032)

ML-135948

ML-136948

Ch.6,0 (0.035)

ML-135960

ML-136960

Ch.7,0 (0.038)

ML-135970

ML-136970

Ch.8,0 (0.038)

ML-135980

ML-136980

Ch.9,0 (0.038)

ML-135990

ML-136990

mediland GmbH • Im Fuchshau 22 • D-73635 Rudersberg (Germany)

Tel.: +49 (0) 7183 30 78 01 0 • Fax: +49 (0) 7183 30 78 01 11 • www.mediland.de • info@mediland.de

"EINE SICHERE PARTNERSCHAFT
FÜR IHRE GESUNDHEIT"

mediland
products & more




Uretersplint


- Polyurethan (PUR)
- super x-ray
- Führungsdraht beschichtet flexible Spitze
- mit Pusher und Fixierklemmen
- steril

Ureteral splint


- Polyurethane (PUR)
- super x-ray
- guide-wire coated flexible tip
- with pusher an fastening clips
- sterile




- Modell ML-133
- Standard
- zylindrische Spitze
- Ch.4,8 / 6,0 / 7,0 / 8,0
- cm 22 / 24 / 26 / 28



- Modell ML-135
- Standard
- zentrale Spitze
- Ch.4,8 / 6,0 / 7,0 / 8,0
- cm 22 / 24 / 26 / 28



- Modell ML-134
- Steuerbar
- zylindrische Spitze
- Ch.4,8 / 6,0 / 7,0 / 8,0
- cm 22 / 24 / 26 / 28



- Modell ML-136
- Steuerbar
- konische Spitze
- Ch.4,8 / 6,0 / 7,0 / 8,0
- cm 22 / 24 / 26 / 28



- Modell ML-138
- Tumor
- zylindrische Spitze
- Ch.4,8 / 6,0 / 7,0 / 8,0
- cm 22 / 24 / 26 / 28



- Modell ML-139
- Tumor
- zentrale Spitze
- Ch.4,8 / 6,0 / 7,0 / 8,0
- cm 22 / 24 / 26 / 28