

"EINE SICHERE PARTNERSCHAFT
FÜR IHRE GESUNDHEIT"

mediland
products & more


products & more
mediland



Silikon Hämaturie Spülkatheter

- 3-Wege Spülkatheter
- Länge 40cm
- zylindrische Spitze
- transparent
- x-ray
- steril
- VE 5 Stück
- UK 40 Stück

ML-1156018 Ch.18 60ml
ML-1156020 Ch.20 80ml
ML-1156022 Ch.22 80ml
ML-1156024 Ch.24 80ml



Silicone haematuria irrigation catheter

- 2 way irrigation catheter
- length 40cm
- cylindrical tip
- transparent
- x-ray
- sterile
- unit 5 pieces
- box 40 pieces

mediland GmbH • Gewerbepark Rudersberg • Ziegeleistrasse 4 • D-73635 Rudersberg (Germany)
Tel.: +49 (0) 7183 30 78 01 0 • Fax: +49 (0) 7183 30 78 01 11 • www.mediland.de • info@mediland.de

"EINE SICHERE PARTNERSCHAFT
FÜR IHRE GESUNDHEIT"

mediland
products & more


products & more
mediland



Silikon Hämaturie Spülkatheter

- 3-Wege Spülkatheter
- Länge 40cm
- Mercier Spitze
- transparent
- x-ray
- steril
- VE 5 Stück
- UK 40 Stück

ML-1156118 Ch.18 60ml
ML-1156120 Ch.20 80ml
ML-1156122 Ch.22 80ml
ML-1156124 Ch.24 80ml



Silicone haematuria irrigation catheter

- 2 way irrigation catheter
- length 40cm
- Mercier tip
- transparent
- x-ray
- sterile
- unit 5 pieces
- box 40 pieces

mediland GmbH • Gewerbepark Rudersberg • Ziegeleistrasse 4 • D-73635 Rudersberg (Germany)
Tel.: +49 (0) 7183 30 78 01 0 • Fax: +49 (0) 7183 30 78 01 11 • www.mediland.de • info@mediland.de

"EINE SICHERE PARTNERSCHAFT
FÜR IHRE GESUNDHEIT"

mediland
products & more

products & more

mediland



Silikon Hämaturie Spülkatheter

- 3-Wege Spülkatheter
- Länge 40cm
- Dufour Spitze
- transparent
- x-ray
- steril
- VE 5 Stück
- UK 40 Stück

ML-1156218 Ch.18 60ml
ML-1156220 Ch.20 80ml
ML-1156222 Ch.22 80ml
ML-1156224 Ch.24 80ml



Silicone haematuria irrigation catheter

- 2 way irrigation catheter
- length 40cm
- Dufour tip
- transparent
- x-ray
- sterile
- unit 5 pieces
- box 40 pieces

mediland GmbH • Gewerbepark Rudersberg • Ziegeleistrasse 4 • D-73635 Rudersberg (Germany)
Tel.: +49 (0) 7183 30 78 01 0 • Fax: +49 (0) 7183 30 78 01 11 • www.mediland.de • info@mediland.de

"EINE SICHERE PARTNERSCHAFT
FÜR IHRE GESUNDHEIT"

mediland
products & more

products & more
mediland



Silikon Hämaturie Spülkatheter


- 3-Wege Spülkatheter
- Länge 40cm
- Couvelaire Spitze
- transparent
- x-ray
- steril
- VE 5 Stück
- UK 40 Stück

ML-1156318 Ch.18 60ml

ML-1156320 Ch.20 80ml

ML-1156322 Ch.22 80ml

ML-1156324 Ch.24 80ml



Silicone haematuria irrigation catheter

- 2 way irrigation catheter
- length 40cm
- Couvelaire tip
- transparent
- x-ray
- sterile
- unit 5 pieces
- box 40 pieces

mediland GmbH • Gewerbepark Rudersberg • Ziegeleistrasse 4 • D-73635 Rudersberg (Germany)
Tel.: +49 (0) 7183 30 78 01 0 • Fax: +49 (0) 7183 30 78 01 11 • www.mediland.de • info@mediland.de