

"EINE SICHERE PARTNERSCHAFT
FÜR IHRE GESUNDHEIT"

mediland
products & more



Anästhesie Atemmaske

Art.Nr.	Größe size	Gewicht weight	Totraum dead space	Farbe colors	Konnektor connector	VE item per box
KM-200	Extra Large - Gr.7	59,5	310 ml	transparent	22 F	50
KM-201	Large - Gr.6	55,5	250 ml	red	22 F	50
KM-202	Medium - Gr.5	51,0	220 ml	yellow	22 F	50
KM-203	Small - Gr.4	37,5	180 ml	green	22 F	50
KM-204	Pediatric - Gr.3	27,0	85 ml	white	22 F	50
KM-205	Infant - Gr.2	17,5	45 ml	-	15 M	50
KM-206	Neonatal - Gr.1	14,5	25 ml	-	15 M	50
KM-207	Pre-Neonatal - Gr.0	11,0	20 ml	-	15 M	50

mediland GmbH • Im Fuchshau 22 • D-73635 Rudersberg (Germany)
Tel.: +49 (0) 7183 30 78 01 0 • Fax: +49 (0) 7183 30 78 01 11 • www.mediland.de • info@mediland.de