


"EINE SICHERE PARTNERSCHAFT  
FÜR IHRE GESUNDHEIT"

**mediland**  
products & more



### Latexschlauch

- dauer temperatur beständig bis +120°C
- auf Rollen zu je 10m
- unsteril
- VE 1 Rolle
- UK 1 Rolle

### Stomach tube

- constantly resistant to temperatures to +120°C
- on reels of 10m length
- non-sterile
- unit 1 reel
- box 1 reel

	<b>I.D.</b>	<b>O.D.</b>
<b>ML-885001</b>	2,0mm	4,0mm
<b>ML-885002</b>	3,0mm	5,0mm
<b>ML-885003</b>	4,0mm	6,0mm
<b>ML-885004</b>	4,0mm	7,0mm
<b>ML-885005</b>	5,0mm	8,0mm
<b>ML-885006</b>	5,0mm	9,0mm
<b>ML-885007</b>	6,0mm	9,0mm
<b>ML-885008</b>	6,0mm	10,0mm
<b>ML-885009</b>	7,0mm	10,0mm
<b>ML-885010</b>	8,0mm	12,0mm
<b>ML-885011</b>	9,0mm	13,0mm
<b>ML-885012</b>	10,0mm	14,0mm

mediland GmbH • Im Fuchshau 22 • D-73635 Rudersberg (Germany)

Tel.: +49 (0) 7183 30 78 01 0 • Fax: +49 (0) 7183 30 78 01 11 • [www.mediland.de](http://www.mediland.de) • [info@mediland.de](mailto:info@mediland.de)