


"EINE SICHERE PARTNERSCHAFT
FÜR IHRE GESUNDHEIT"

mediland
products & more



Bronchial Doppellumige Tracheostomiekannüle

- Großvolumige Niederdruck-Cuffs
- Wolfram-Ringe zur röntgenologischen Lagekontrolle
- patentierter strikturfreier Übergang von den Konnektoren zum Doppellumenschlauch für die problemlose Passage eines Bronchoskopes
- effektives Verhältnis vom inneren zum äußeren Querschnitt
- steril
- VE 1 Stück
- UK 1 Stück

Zubehör zur „DLTK“:

- Absaugkatheter mit Saugregulierung
- Winkelanschlußstücke
- Konnektorschlüssel

Endo bronchial double lumen tracheostomy cannula

- high volume low pressure cuffs
- stricture free transition between single lumen and double lumen
- patented for effective passage of a bronchoscope
- x-ray markers of tungsten
- effective relation between inside and outside cross-section
- -sterile
- unit 1 piece
- box 1 pieces

Accessoires for the „DLTK“:

- suction catheter with finger valve
- connecting pieces (angled)
- disconnecting device (wrench)

links / left	O.D.	intratracheale Länge/intracheal length	rechts / right
ML-11643585	Ch.35	85mm	ML-11653585
ML-11643775	Ch.37	75mm	ML-11653775
ML-11643785	Ch.37	85mm	ML-11653785
ML-11643795	Ch.37	95mm	ML-11653795
ML-11643975	Ch.39	75mm	ML-11653975
ML-11643985	Ch.39	85mm	ML-11653985
ML-11643995	Ch.39	95mm	ML-11653995
ML-11644195	Ch.41	95mm	ML-11654195

mediland GmbH • Im Fuchshau 22 • D-73635 Rudersberg (Germany)

Tel.: +49 (0) 7183 30 78 01 0 • Fax: +49 (0) 7183 30 78 01 11 • www.mediland.de • info@mediland.de

